



DVH Analytics



Application Views

- Bin-wise statistics calculated (median, mean, quartiles)
- Interactive plot: pan, zoom, mouse-over info
- Single review DVH calculated on the fly from a temporary folder (not from database).  
Use to compare a current plan to historical data

DVHS

MRN to review


ANON99785

ROI to review

BrainStem

Rx Dose (Gy):

52.08





# DVH Analytics

Application Views





# DVH Analytics

## Application Views

### ROI Viewer

- Provides a simple axial slice of ROI data (no CT image)
- MRNs restricted to query, but all ROIs for the selected patient are available
- Currently, only 5 ROIs at a time
- First ROI defaults to external/skin/body
- Plan to add distance measurement tools